

Description

Features:

1. Ultrasonic heating function increase blood circulation, slow down the cell excited point, relief pain.
2. Makes the subcutaneous tissue movement, activation of cells, reduce puffiness of the face, smooth tissue congestion, inhibition of dark mark formation, removal.
3. Improvement of microcirculation, inhibiting facial varicose, improve dry skin and skin aging lines, detox.
4. Loosen the nerve muscle, resist muscular nerve.
5. Dissolve hypodermic fat, accelerate absorption, heals acne skin.
6. Light weight and compact design.

Main Functions:

1. Significant effect on removes wrinkles, pimples and spots
2. Skin care to clear, lighten, smooth the look of skin
3. Enhances blood circulation
4. Deep cleanse and moisturize
5. Enhance skin texture, anti-wrinkle, helps fight fatty areas, tightens and firms

In only 10-15 minutes daily, you'll see and experience dynamic changes in your face: your skin will appear more youthful, smooth and healthy with this intense form of care. Ultrasound Beauty Massagers reach deeply into the skin. Although invisible to the eye and undetectable by our skin, Ultrasound Beauty Devices radiate micro-vibrations, which are often used to help stimulate the body's recovery from illness and injury. Waves of ultrasound can also reach fats, built up in areas such as lower abdomen, hips, thighs, calves or double chin.

Specification:

Color:White

Voltage: AC 220V or 110V UK EU US AU Plug

Power:10W

Adjustable timer, adjustable intensity

Treatment area: Facial, body

Product size: 40.5x27.5x11cm

About the Power Selection(Level Set): Usually Level 1-3 is for eye,Level 4-6 is for face and Level 6 is for body.For sensitive skin,please select the low level.

About the Ultrasonic Wave Selection(Type Set): The light keeps on is usually for importing skin care products; the light flash slowly is for eye care; the light flash quickly is for exporting toxin.

Please don't export toxin around eye area.

Package Included:

1X Main Ultrasound Unit

1X Face Ultrasound Probe

1X Eye Ultrasound Probe

Note:Intrustction only available in Chinese.